

Lab  
Power

# IPA DISINFECTANT WIPES



  
**westlab**  
Enriching the World of Science

# ISOPROPYL ALCOHOL DISINFECTANT WIPES

664-231

## DESCRIPTION

The first 'Australian Made' wipe by Westlab. A hospital grade wipe manufactured here in Australia. Hospital Grade.

Available for export with current production capacity of 50,000 canisters per week.

Wall brackets to suit, white powder coated bracket that mounts to the wall.

Made in Australia

## SPECIFICATIONS

Isopropyl Alcohol	70%
Wipes per canister	150 wipes
Wipe dimension	200mm x 150mm
Wipe material	Thermalbond polypropylene
Wipe thickness	17gsm
Origin	Made in Australia
Canister dimension	100mm x 173mm
Canisters per carton	12 canisters
Canisters per pallet	648 canisters

## ACCESORIES

Wall Bracket for LabPower  
IPA Disinfectant Wipes

Code: 664-232

