

# Finger Pulse Oximeter

## XY301

Easy to use, light, and portable oximeter.  
For hospital, homecare, sport and etc.

### Product Features

- Automatic shutdown after 8s of no signal
- Three-axis accelerometer
- Pulse wave display
- Pulse column display

### Product Specifications

Power consumption	≤30mA
Humidity	≤80%
Boot method	Push button to boot
Power supply	2 x 1.5V AAA alkaline batteries
Display mode	OLED
Switch of display direction	4 directions
QTN	40PCS/CTN
G.W./N.W.	3.8/1.2kgs
Size	45cm x 26cm x 15cm
ARTG number	380874

For additional information please contact Westlab on 1800 358 101 or go to [www.westlab.com.au](http://www.westlab.com.au).

Small volume,  
low power  
consumption,  
convenient size.

