

ISOLATION GOWN


Level 2, Non-sterile



TIE ON THE BACK NECK



LONG SLEEVE, KNITTED CUFFS



TIE AROUND THE WAIST



ISOLATION GOWN LEVEL 2 NON-STERILE

KEY FEATURES

Material: PP with PE Clothing, PP non-woven

Weight: PP with PE Coating – 18 grams PP + 12 grams PE +/- 2 grams | PP non woven – 25 grams +/- 2 grams

Size: 120 x 140 cm

Colour: PP with PE Coating – Blue, PP non woven – Yellow

Standard: EN14126, Class 2/Class 3/Class 4

Application: Single Use, General Medical, Labs, Patient Care & Critical Care Units

Package: 10 pcs/bag, 10 bags/ctn, 365 x 280 x 450mm

Features: Knitted Cuff, Tie around waist, tie around back of neck

SIZING

Size	S	M	L	XL	XXL	XXXL	Tolerance
A	110	115	127	135	140	145	+2
B	140	145	150	150	155	160	+2
C	19	20	21	22	23	24	±1
D	61	63	65	68	70	73	+1
E	18	18	18	18	19	20	+1
F	28	29	29	32	32	33	+1
I	59	59	73	73	77	80	+1
L	100	105	117	125	130	135	±1
M	55	58	60	60	64	66	+1
N	3'9	3'9	3'9	3'9	3'9	3'9	+1
O	3'6	3'6	3'6	3'6	3'6	3'6	+1

ORDERING WITH WESTLAB IS EASY...

ONLINE

www.westlab.com.au

Order online with our easy Express
Order form for fast, secure ordering.



PHONE

1800 385 101

Speak to one of our
friendly, efficient staff
from 8.00am to 5.00pm Monday to
Friday.



EMAIL

sales@westlab.com.au

24 hours, 7 days a week convenience.



FAX

1300 725 903

Order via fax and we will confirm your
order and delivery details at the end
of the day.

